medi foot sport

For active sports

- shock absorption
- high endurance
- breathable

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medi foot sport profi

For professional sports

- extra shock absorption
- high endurance
- breathable

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medi foot comfort

For sensitive feet

- gentle arch support
- makes every shoes comfortable

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medi foot light

For women's fancy shoes

- ultra thin and lightweight
- for every heel height

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medi foot travel

For high performance

- rigid arch support
- protection of feet in cold exposure

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medi foot soft

For problem feet

- gentle support
- correction of foot defects
- keep feet at a comfortable temperature
- uniformly distributed load

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Orthopaedic insoles

Dear Patient,

Congratulations! You have received a pair of medi foot orthopaedic insoles.

medi foot product line is a result of collaboration between, technologists, biomechanics and designers.

In medical statistics, 8 out of 10 people will develop or already have flat feet. This may lead to inadequate shock absorption, joint and vertebral deformation. Early signs of flat feet may be Heavy Leg Syndrome, feet/shank pain, swollen feet and limitation of movement.

medi foot orthopaedic insoles are innovative functional foot orthoses. When used correctly they will provide many years of comfort and support. Wide range of insoles can be used not only as a treatment of flat feet but to cushion foot impact for a variety of activities. Properly prescribed insoles provide arch support, optimize biomechanics of movement and off-load your feet (e.g. when you do sport activity).

medi foot insoles are made of ecological materials using state-of-the-art technologies and meet all of the requirements of the European standards. In accordance with regulations, our products are subject to the conformity.

This will ensure high medical efficacy and comfort of our products.

Indications

- treatment of flat foot (I-II stage)
- rehabilitation after foot injury
- foot weakness due to functional insufficiency
- extra foot protection during sports activity, pregnancy, diabetes

Contraindications / Side effects

Individual intolerance of materials

Mode of action

- transverse and longitudinal arch support
- extra shock absorption and reduction of pain

Wearing instruction

Insoles should be used with regular closed footwear.

It is normal for the insoles to feel a bit odd at first (till 2 weeks). However, if at any time during the break-in period you experience pain that lasts for more than a day in your ankles, knees, hips, or back you should stop wearing the insoles and contact your doctor.

Cleaning instructions

Clean with a damp cloth using a mild detergent. Never use ointment, lotion or oil as detergent. Allow insoles to air dry (20-25 °C). Do not air dry near a heat source. Do not use solvents for remove stains. Do not dry-clean.

Storage instructions

Store in a cool dry place and protect from direct sunlight. 6 month limited warranty.

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